

Product datasheet for MR201972

Sri (NM_001080974) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sri (NM_001080974) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Sri
Synonyms: 2210417O06Rik; 2900070H08Rik; Sor
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR201972 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCGTATCCCGGGCACCCCTGGCGGGCGGAGGGTACTACCCAGGCGGGTATGGTGGGGCTCCCGGAG
 GGCCTGCGTTCCTGGACAACTCAGGATCCGCTATATGGTTACTTTGCTGCTGTGGCTGGACAGGACGG
 ACAAAATTGATGCTGACGAATTGCAGAGATGCCTAACACAGTCTGGCATTGCGGGAGGATACAAGCCTTTT
 AACCTGGAGACTTGTGCGCTTATGGTTCAATGTTGGATAGAGATATGTCGGCACGATGGGGTTCAATG
 AATTCAAAGAGCTCTGGGCTGTGCTGAATGGCTGGAGACAACACTTCATCAGTTTCGATAGCGATAGGAG
 TGGCACCGTGGATCCCCAGGAGCTACAGAAGGCCCTGACGACTATGGGATTCAGGTTGAGCCCCAACT
 GTGAATTCAGTTGCAAAGCGATACAGCACCAGCGGAAGATCACGTTTGATGACTACATTGCCTGCTGTG
 TCAAACAGGGCTCTCACAGATAGCTTTCGACAGCGGGATTCGGGTCAACAAGGAGTTGTGAATTTCTC
 GTATGATGATTTTCATCCAGTGTGTCATGACCGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR201972 protein sequence
 Red=Cloning site Green=Tags(s)

MAYPGHPGAGGGYYPGGYGGAPGGPAFPGQTQDPLYGYFAAVAGQDQIDADELQRCLTQSGIAGGYKPF
 NLETCLMVSMLDRDMSGTMGFNEFKELWAVLNGWRQHFISFSDRSRGTVDPQELQKALTTMGFRLSPQT
 VNSVAKRYSTSGKITFDDYIACCVKLRALTDSFRRRDSGQGVVNF SYDDFIQCVMTV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001080974

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080974.2](#), [NP_001074443.1](#)

RefSeq Size: 2399 bp

RefSeq ORF: 597 bp

Locus ID: 109552


UniProt ID: [Q6P069](#)

Cytogenetics: 5 3.38 cM

MW: 21.6 kDa

Gene Summary: Calcium-binding protein that modulates excitation-contraction coupling in the heart. Contributes to calcium homeostasis in the heart sarcoplasmic reticulum. Modulates the activity of RYR2 calcium channels.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR201972