

Product datasheet for **MR201989L3V**

Commd9 (NM_029635) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Commd9 (NM_029635) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Commd9
Synonyms:	1810029F08Rik; AA409293
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_029635
ORF Size:	594 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201989).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_029635.3 , NP_083911.1
RefSeq Size:	1175 bp
RefSeq ORF:	597 bp
Locus ID:	76501
UniProt ID:	Q8K2Q0
Cytogenetics:	2 E2



[View online »](#)

Gene Summary:

May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes. May down-regulate activation of NF-kappa-B. Modulates Na(+) transport in epithelial cells by regulation of apical cell surface expression of amiloride-sensitive sodium channel (ENaC) subunits.
[UniProtKB/Swiss-Prot Function]