

Product datasheet for **MR201992L4V**

Prdx1 (NM_011034) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Prdx1 (NM_011034) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Prdx1
Synonyms:	MSP23; NkefA; OSF-3; OSF3; PAG; Paga; Prdxl; prx1; Prxl; TdpX2; TDX2; TPxA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_011034
ORF Size:	600 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201992).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011034.2
RefSeq Size:	1417 bp
RefSeq ORF:	600 bp
Locus ID:	18477
UniProt ID:	P35700
Cytogenetics:	4 53.28 cM



[View online »](#)

Gene Summary:

Thiol-specific peroxidase that catalyzes the reduction of hydrogen peroxide and organic hydroperoxides to water and alcohols, respectively. Plays a role in cell protection against oxidative stress by detoxifying peroxides and as sensor of hydrogen peroxide-mediated signaling events. Might participate in the signaling cascades of growth factors and tumor necrosis factor-alpha by regulating the intracellular concentrations of H₂O₂ (By similarity). Reduces an intramolecular disulfide bond in GDPD5 that gates the ability to GDPD5 to drive postmitotic motor neuron differentiation (PubMed:19766572).[UniProtKB/Swiss-Prot Function]