

Product datasheet for **MR201998L4V**

Zmat2 (NM_025594) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zmat2 (NM_025594) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Zmat2
Synonyms:	2610510D14Rik; 2900082I05Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_025594
ORF Size:	597 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR201998).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025594.3 , NP_079870.1
RefSeq Size:	2257 bp
RefSeq ORF:	600 bp
Locus ID:	66492
UniProt ID:	Q9CPW7
Cytogenetics:	18 B2
Gene Summary:	Involved in pre-mRNA splicing as a component of the spliceosome.[UniProtKB/Swiss-Prot Function]



[View online »](#)