

Product datasheet for **MR202000L3V**

Ung (BC004037) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ung (BC004037) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ung
Synonyms:	UNG1, UNG2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC004037
ORF Size:	599 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202000).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC004037
RefSeq Size:	1737 bp
RefSeq ORF:	599 bp
Locus ID:	22256
Cytogenetics:	5 F
Gene Summary:	Excises uracil residues from the DNA which can arise as a result of misincorporation of dUMP residues by DNA polymerase or due to deamination of cytosine.[UniProtKB/Swiss-Prot Function]



[View online »](#)