

Product datasheet for **MR202005L4V**

Pemt (NM_008819) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pemt (NM_008819) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pemt
Synonyms:	PEAMT; Pempt; Pempt2; PEMT2; PLMT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_008819
ORF Size:	600 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202005).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_008819.2
RefSeq Size:	1023 bp
RefSeq ORF:	600 bp
Locus ID:	18618
UniProt ID:	Q61907
Cytogenetics:	11 37.81 cM



[View online »](#)

Gene Summary:

Catalyzes the three sequential steps of the methylation pathway of phosphatidylcholine biosynthesis, the SAM-dependent methylation of phosphatidylethanolamine (PE) to phosphatidylmonomethylethanolamine (PMME), PMME to phosphatidylmethylethanolamine (PDME), and PDME to phosphatidylcholine (PC). [UniProtKB/Swiss-Prot Function]