

Product datasheet for **MR202016L3V**

Rpp25 (NM_133982) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Symbol: Rpp25

Synonyms: A1851155

Mammalian Cell Selection: Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK

ACCN: NM_133982

ORF Size: 597 bp

ORF Nucleotide Sequence: The ORF insert of this clone is exactly the same as(MR202016).

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

RefSeq: [NM_133982.1](#), [NP_598743.1](#)

RefSeq Size: 1361 bp

RefSeq ORF: 600 bp

Locus ID: 102614

UniProt ID: [Q91WE3](#)

Cytogenetics: 9 B



Gene Summary:

Component of ribonuclease P, a ribonucleoprotein complex that generates mature tRNA molecules by cleaving their 5'-ends. Also a component of the MRP ribonuclease complex, which cleaves pre-rRNA sequences.[UniProtKB/Swiss-Prot Function]