

Product datasheet for **MR202035L3V**

Commd7 (NM_133850) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Commd7 (NM_133850) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Commd7
Synonyms:	2310010I22Rik; AW060956; mU3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_133850
ORF Size:	600 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202035).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133850.2 , NP_598611.2
RefSeq Size:	3213 bp
RefSeq ORF:	603 bp
Locus ID:	99311
UniProt ID:	Q8BG94
Cytogenetics:	2 H1
Gene Summary:	May modulate activity of cullin-RING E3 ubiquitin ligase (CRL) complexes. Associates with the NF-kappa-B complex and suppresses its transcriptional activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)