

Product datasheet for **MR202074A1V**

Mouse Pmf1 (NM_025928) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Pmf1 (NM_025928) AAV Particle
Tag: Myc-DDK
Symbol: Pmf1
Synonyms: 2600009M07Rik; AL033286; AW060657
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202074 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTGAGGTGAGCCGAGATAGCGAGGCTGCGGAAAGGGGGCCTGAGGGCTCCTCTCCGAAGCTGTGC
CAGGGGACGCGACCATCCCCAGGGTGAACTCCTGGACGCCATAGTAGACACTTTCCTCCAGAAGCTAGT
CGCCGACAGGAGCTACGAGAGGTTACACACCTGCTACAAACACTTCCACCAGTTGAACCTGAGGTGACG
CAGAGGATCTATGACAAGTTTGTGGCTCAGTTGCAGACATCCATCCGCGAGGAAATCTCAGAAATCAAAG
AGGAGGGGAACCTAGAAGCTGTCTGAACCTCCCTGGATAAGATCATAGAAGAAGGCAGAGAGCGCGGAGA
GCCAGCCTGGCGACCCAGTGGAATCCCAGAGAAAGACCTGTGTAGTGTGCATGGCACCCCTACTTCTGAAG
CAACAGGATACCTGTGTGCATCAAGTACGGAACAGGAAGCCAAGAACCAGGAAGTGGCCGACGCTGTCC
TGGCCGGGCGCAGGCAGGTGGAGGAGCTGCAGCAGCAGGTTCCGGGCCCTCCAGCAAAATGGCAGGCTCT
ACACAGAGAGCAGAGGGAGCTGCTGTCTAGTGTGAGGGCGCCTGAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202074 protein sequence
Red=Cloning site Green=Tags(s)

MAEVSRDSEAAERGPEGSSPEAVPGDATIPRVKLLDAIVDTFLQKL VADRSYERFTTCYKHFHQLNPEVT
QRIYDKFVAQLQTSIREEISEIKEEGNLEAVLNSLDKIIIEGRERGEPAWRPSGIPEKDLCSVMAPYFLK
QQDTLCHQVRKQEAQNQELADAVLAGRRQVEELQQQVRALQQTWQALHREQRELLSVLRAPE

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse



[View online »](#)

Serotype:	AAV-2
ACCN:	NM_025928
ORF Size:	606 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_025928.1</u>
RefSeq Size:	992 bp
RefSeq ORF:	609 bp
Locus ID:	67037
UniProt ID:	<u>Q9CPV5</u>
Cytogenetics:	3 F1
MW:	23.1 kDa