

Product datasheet for MR202074L4

Pmf1 (NM_025928) Mouse Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pmf1 (NM_025928) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Pmf1
Synonyms:	2600009M07Rik; AL033286; AW060657
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202074).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_025928
ORF Size:	609 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_025928.1
RefSeq Size:	992 bp
RefSeq ORF:	609 bp
Locus ID:	67037
UniProt ID:	Q9CPV5
Cytogenetics:	3 F1
Gene Summary:	Part of the MIS12 complex which is required for normal chromosome alignment and segregation and for kinetochore formation during mitosis (By similarity). May act as a cotranscription partner of NFE2L2 involved in regulation of polyamine-induced transcription of SSAT.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR202074L4