

Product datasheet for **MR202077L4V**

Ccdc85b (NM_198616) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ccdc85b (NM_198616) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ccdc85b
Synonyms:	AI842788; BB115872
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_198616
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202077).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_198616.1 , NP_941018.1
RefSeq Size:	4401 bp
RefSeq ORF:	609 bp
Locus ID:	240514
UniProt ID:	Q6PDY0
Cytogenetics:	19 A



[View online »](#)

Gene Summary:

Functions as a transcriptional repressor. May inhibit the activity of CTNNB1 in a TP53-dependent manner and thus regulate cell growth. May function in adipocyte differentiation, negatively regulating mitotic clonal expansion (PubMed:15644333, PubMed:22666460). [UniProtKB/Swiss-Prot Function]