

Product datasheet for **MR202095L4V**

Edn1 (NM_010104) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Edn1 (NM_010104) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Edn1
Synonyms:	ET-1; PPET1; preproET
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_010104
ORF Size:	609 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202095).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010104.2
RefSeq Size:	2344 bp
RefSeq ORF:	609 bp
Locus ID:	13614
UniProt ID:	P22387
Cytogenetics:	13 20.82 cM



[View online »](#)

Gene Summary:

This gene encodes a member of the endothelin family of peptides. The encoded preproprotein undergoes proteolytic processing to generate a peptide before secretion by the vascular endothelial cells. The mature peptide has various biological activities such as vasoconstriction, cell proliferation, stimulation of hormone release and modulation of central nervous activity. Mice lacking the encoded protein exhibit neonatal lethality accompanied with numerous craniofacial and cardiovascular defects due to disruption in cranial and cardiac neural crest cell patterning during early embryogenesis. [provided by RefSeq, Feb 2016]