

Product datasheet for **MR202116L3V**

Asf1a (NM_025541) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Asf1a (NM_025541) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Asf1a |
| Synonyms: | 2310079C17Rik; C79052 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_025541 |
| ORF Size: | 615 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR202116). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_025541.3 , NP_079817.1 |
| RefSeq Size: | 2192 bp |
| RefSeq ORF: | 615 bp |
| Locus ID: | 66403 |
| UniProt ID: | Q9CQE6 |
| Cytogenetics: | 10 B3 |



[View online »](#)

Gene Summary:

Histone chaperone that facilitates histone deposition and histone exchange and removal during nucleosome assembly and disassembly. Cooperates with chromatin assembly factor 1 (CAF-1) to promote replication-dependent chromatin assembly and with HIRA to promote replication-independent chromatin assembly. Required for the formation of senescence-associated heterochromatin foci (SAHF) and efficient senescence-associated cell cycle exit. [UniProtKB/Swiss-Prot Function]