

Product datasheet for **MR202131**

Polr3h (NM_030229) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Polr3h (NM_030229) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Polr3h
Synonyms:	5031409G22Rik; RPC8
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202131 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTCTGCTGGTGGAGATGGTGGACACCGTGCGCATACCCCCGTGGCAGTTTGAGCGGAACTCAACG
ATTCCATTGCTGAGGAGCTGAACAAGAAATGGCCAACAAGGTCGTATACAACGTGGGACTCTGCATCTG
TCTGTTTGATATACCAAGCTGGAGGACGCTATGTGTTCCAGGGGATGGAGCATCCCACCAAAGTC
CATTTCCGCTACGTGGTGTTCACCCCTCCTGGATGAGATCCTATTGGGAAGATCAAAGGCTGCAGCC
CAGAGGGTGTGCACGTCTTTAGGGTCTTTGATGACATTCTCATCCCCCTGAGTCACTGCAGCAGCC
TGCCAAGTTTGATGAAGCGGAGCAAGTGTGGGTGTGGGAGTACGAGACGGAGGAAGGAGCGCATGACCTG
TACATGGACACTGGGAGGAGATCCGCTTCAGGGTGGTGGATGAGAGCTTTGTTGGACACATCCCCACAG
GACCCAGTTCCGCCGAGGCTGCCTCTTCAAGCGAGGAGCTGCCAAGAAGGAAGCACCATATACACTTGT
GGGATCTATCAGTGAGCCAGGACTAGGCCTTCTCTCTGGTGGACGAGCAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>MR202131 protein sequence Red=Cloning site Green=Tags(s)
-------------------	--

MFVLVEMVDTVRIPPWQFERKLNDSIAEELNKKLANKVVVNVGLCICLFDITKLEDAYVFPDGDASHTKV
HFRYVVFHPFLDEILIGKIKGCSPEGVHVS LGFFDDILIPPELQQPAKFDEAEQVWWEYETEEGAHDL
YMDTGEEIRFRVDES FVDTSP TGPSSAE AASSSEELPKKEAPYTLVGSISEPGLLLSWWTSN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_030229

ORF Size: 615 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_030229.4](#), [NP_084505.2](#)

RefSeq Size: 2514 bp

RefSeq ORF: 615 bp

Locus ID: 78929


UniProt ID: [Q9D2C6](#)

Cytogenetics: 15 E1

MW: 22.9 kDa

Gene Summary: DNA-dependent RNA polymerase catalyzes the transcription of DNA into RNA using the four ribonucleoside triphosphates as substrates. Specific peripheric component of RNA polymerase III which synthesizes small RNAs, such as 5S rRNA and tRNA. Plays a key role in sensing and limiting infection by intracellular bacteria and DNA viruses. Acts as nuclear and cytosolic DNA sensor involved in innate immune response. Can sense non-self dsDNA that serves as template for transcription into dsRNA. The non-self RNA polymerase III transcripts induce type I interferon and NF- Kappa-B through the RIG-I pathway (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR202131