

Product datasheet for MR202137

Rps5 (NM_009095) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rps5 (NM_009095) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rps5
Synonyms:	AA617411
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202137 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACTGAGTGGGAAGCAGCCACACCAGCGGTGGCAGAGACCCCTGACATCAAGCTCTTTGGGAAATGGA
GCACTGATGACGTGCAGATCAACGATATTTCTCTGCAGGATTACATTGCTGTGAAGGAGAAGTATGCCAA
GTACCTGCCCCACAGTGCCGGACGGTATGCTGCCAAGCGCTTCGCAAAGCACAATGTCCCATCGTGGAG
CGCCTTACTAACTCCATGATGATGCATGGTCGTAACAACGGCAAGAAGCTCATGACTGTGCGAATTGTCA
AGCATGCCTTTGAGATCATCCACCTGCTCACTGGTGAGAACCCCTCTGCAGGTCCTGGTGAATGCTATCAT
CAACAGTGGCCCCGAGAAGACTCAACACGCATTGGGCGGGCCGGTACAGTGAGACGACAGGCTGTGGAT
GTGTCCCCACTGCGTCGAGTGAATCAGGCCATCTGGCTGCTGTGCACAGGGGCTCGTGAGGCTGCTTTCC
GGAACATCAAGACCATCGCCGAGTGCCTTGAGATGAGCTCATTAAATGCTGCCAAGGGCTCCTCCAATTC
CTATGCCATCAAGAAGAAAGATGAAGTGGAGCGTGTGGCCAAGTCTAACC GC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:	>MR202137 protein sequence Red=Cloning site Green=Tags(s)
-------------------	--


MTEWAATPAVAETPDIKLFGKWSTDDVQINDISLQDYIAVKEKYAKYLPHSAGRYAAKRFKAQCPIVE
RLTNSMMMHGRNNGKKLMTVRIVKHAFIHHLLTGENPLQVLVNAIINSGPREDSTRIGRAGTVRRQAVD
VSPLRVNQAIWLLCTGAREAAFRNIKTIAECLADELINAAGSSNSYAIKKKDELERVAKSNR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_009095

ORF Size: 612 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.


RefSeq: [NM_009095.2](#)

RefSeq Size: 768 bp

RefSeq ORF: 615 bp

Locus ID: 20103
 UniProt ID: [P97461](#)
 Cytogenetics: 7 A1
 MW: 22.9 kDa

Product images:



Circular map for MR202137