

Product datasheet for MR202152

Rhox4a (NM_001039688) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rhox4a (NM_001039688) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Rhox4a
Synonyms: Rhox4; Rhox4.1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR202152 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGAGCATCAAACACCAACTACCTACTTCATGAGGGACTTGGCAAGGACAAGGAAAATTTGAATGGTG
 GGAAGACACAGGCAGTCTTACCACTGGATGGAGAGGGAAGAAATGAGGGAGAGAGTACTGGGCCAGTC
 CGGAGCCCGCAGCAGTGAAGGGGACAAAGCAGAAGAATTAAGTGGAGAAGGTGGGCTGCTGCTGGTGAT
 GCAGACCTCATGGATAACAGCAACCAAGAGGACCAGGACACCAGTGGCAGTGCCAGGAGGAGGAGAAGC
 TGCCAGAGGAGCCAGTTCTCAAGGATGCTGTGGTCATAGACAAAGTGCAGCCTATTCCAGTCTGGTATC
 TGGTGTGCGGCCAAGTCAGTGTGGGTACAGCAGCGTAGCTTACACTACAATTTCCAATGGTGGCAGCTG
 CAGGAGCTGGAGCGCATTTCCAGCAGAATCACTTCATCCGTGCAGAGGAAAGAAGACATCTGGCAAGAT
 GGATAGGTGTGAGTGAAGCCAGAGTTAAGAGATGGTTTAAAGAAGAGGAGAGACACTTCAGAAGAGGACA
 AAGTCAGTTAGGAATGAATGATGACGCCTCTGTGGGTCCCACTCTACCTTTCTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202152 protein sequence
 Red=Cloning site Green=Tags(s)

MEHQNTNYLLHEGLGKDKENLNGGKTQAVLPLDGEGRNEGESVLGQSGAAAVEGDKAEELSGEGGPAAGD
 ADLMDNSNQEDQDTSGSAQEELPEEPVLKDAVVIDKVQPIPVLVSGVRPKSVVWVQQRSLHYNFQWWQL
 QELERIFQQNHFIRAEERRHLARWIGVSEARVKRWFKKRREHFRRGQSQLGMNDDASVGSHTFL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001039688

ORF Size: 618 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:**
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001039688.3](#), [NP_001034777.1](#)

RefSeq Size: 891 bp


RefSeq ORF: 618 bp

Locus ID: 664609

Cytogenetics: X A3.3

MW: 23.1 kDa

Product images:



Circular map for MR202152