

## Product datasheet for **MR202160L3V**

### Mad2l1 (NM\_019499) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Mad2l1 (NM_019499) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mad2l1
Synonyms:	AA673185; MAD2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019499
ORF Size:	615 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202160).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019499.4</a>
RefSeq Size:	1725 bp
RefSeq ORF:	618 bp
Locus ID:	56150
UniProt ID:	<a href="#">Q9Z1B5</a>
Cytogenetics:	6 30.56 cM



[View online »](#)

**Gene Summary:**

Component of the spindle-assembly checkpoint that prevents the onset of anaphase until all chromosomes are properly aligned at the metaphase plate. Required for the execution of the mitotic checkpoint which monitors the process of kinetochore-spindle attachment and inhibits the activity of the anaphase promoting complex by sequestering CDC20 until all chromosomes are aligned at the metaphase plate (By similarity).[UniProtKB/Swiss-Prot Function]