

Product datasheet for MR202169

Ppp1r2 (BC069886) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ppp1r2 (BC069886) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ppp1r2
Synonyms: 0610025N14Rik; 2310007G06Rik; 4930440J04Rik; 5430408E15Rik; D16Ertd248e; IPP-2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR202169 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGGCTCAACGGCCTCGCACCGGCCATCAAAGGAATCCTGAAGAACAAGACCTCTGCGGCTCCC
CCCCGGTGGTGCCTCGGCCGAACAGCCCCGCCGATCGTGAAGAGGAAGTCTGAGTAAGAAGTCTCAGAA
GTGGGATGAAATGAATATCCTGGCAACATATCATCCCGCTGACAAAGACTATGGCTTGATGAAAATAGAC
GAGCCAAACACCCATATCACAATATGATAGGTGACGATGAAGACGCCTATAGTGATTCAGAAGGCAATG
AAGTCATGACTCCAGATATCTTAGCTAAGAAATTAGCTGCTGCTGAAGGCTCAGAGCCAAAGTATCGGAC
TCGGGAGCAAGAAAGTAGTGGAGAGGAAGATAATGACCTCTCCCTGAGGAGCGAGAAAAGAAGCGACAG
TTTGAAATGAAAAGGAAGCTTATTACAATGAAGGACTCAATATTAAGTTAGCCAGACAGTTAATTTCAA
AAGACCTACATGATGATGATGAAGATGAAGAAATGGCAGAGACTGCAGATGGAGACAGCATGAATGTTGA
AGAGTCAAGTCAAGGATCTACTACAAGTGACCATCTGCAACACAAATCACAAAGTTCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202169 protein sequence
Red=Cloning site Green=Tags(s)


MAASTASHRPIKILKNKTSASPPVPSAEQPRPIVEEELSKKSQKWDEMNILATYHPADKDYGLMKID
EPNTPYHNMIGDDEDAYSDEGNEVMTDILAKKLAEEGSEPKYRTREQESSGEEDNDLSPEEREKKRQ
FEMKRKLHYNELNIKLARQLISKDLHDDDEDEEMAETADGDSMNVEESSQGSTTSDHLQHKSQSS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC069886

ORF Size: 618 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC069886.1](#)

RefSeq Size: 1490 bp

RefSeq ORF: 620 bp


Locus ID: 66849

Cytogenetics: 16 21.41 cM

MW: 23.1 kDa

Gene Summary: Inhibitor of protein-phosphatase 1.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR202169