

Product datasheet for **MR202191L3V**

Sike1 (NM_025679) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sike1 (NM_025679) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sike1
Synonyms:	2810005O12Rik; 5730470L24Rik; AI450236; AI839862; Sike; Sikeb
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_025679
ORF Size:	624 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202191).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_025679.2
RefSeq Size:	2798 bp
RefSeq ORF:	624 bp
Locus ID:	66641
UniProt ID:	Q9CPR7
Cytogenetics:	3 F2.2



[View online »](#)

Gene Summary:

Physiological suppressor of IKK-epsilon and TBK1 that plays an inhibitory role in virus- and TLR3-triggered IRF3. Inhibits TLR3-mediated activation of interferon-stimulated response elements (ISRE) and the IFN-beta promoter. May act by disrupting the interactions of IKBKE or TBK1 with TICAM1/TRIF, IRF3 and DDX58/RIG-I. Does not inhibit NF-kappa-B activation pathways (By similarity).[UniProtKB/Swiss-Prot Function]