

Product datasheet for **MR202213L4V**

Fn3k (NM_001038699) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fn3k (NM_001038699) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fn3k
Synonyms:	2310074G21Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001038699
ORF Size:	624 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202213).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001038699.2 , NP_001033788.1
RefSeq Size:	3882 bp
RefSeq ORF:	627 bp
Locus ID:	63828
Cytogenetics:	11 E2
Gene Summary:	May initiate a process leading to the deglycation of fructoselysine and of glycosylated proteins. May play a role in the phosphorylation of 1-deoxy-1-morpholinofructose (DMF), fructoselysine, fructoseglycine, fructose and glycosylated lysozyme (By similarity). [UniProtKB/Swiss-Prot Function]



[View online »](#)