

Product datasheet for **MR202215**

Xlr (NM_011725) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Xlr (NM_011725) Mouse Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Xlr
 Mammalian Cell Selection: Neomycin
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 ORF Nucleotide Sequence: >MR202215 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAAACTGGGACTTGTCAGTGATGAAATGCAAGATGGGAATGCTCCAGAATTGGACGTTATTGAAG
 AACATAATCCAGTAACCTCGTGATGATGAGAATGCAATCCTGAAGAAGTAGTTGGAGATACACGATCTCC
 AGTACAAAATATTCTGGGAAATTTGAAGGTGACATTAAACAGAGGCTTCACATAAGAGAAAACGCATG
 GAACTTATATCAAAGATTCTTTCAAAGACAGCAACGTGAAATTAGAACAACCTTTGGAAAACGAACAAAC
 AAGAGAGAAAGAAGATTAACAACAAATTTGTGAGCAGTATATACTACATTTTCAAGATTTGATATGGA
 TGTCCAGAAATTCAATGAAGAACAAGAAAAATCAGTGAACAATTATCAAAAAGAACAACAAGCATTAAAA
 CTGTCCAAATGTAGTCAGAGCCAGACCCTGGAAGCAATTAAGACATGCATGAAACTACATGGAGGGTC
 TGATGAACCTGGAGACCAACAATAACAATATGCTTTTGTAGTGGTGAAGTGAAGAAAGAAATGTC
 AGTGTTTAAAAAGACCTCATGAAGCACTCTGAAGTACTCTTCTAGTTTCCCATCTTCAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202215 protein sequence
 Red=Cloning site Green=Tags(s)

MENWDLSSDEMQDGNAPELDVIEEHNPTVRDDENANPEEVVGDTRSPVQNILGKFEGDINKRLHIKRKRM
 ETYIKDSFKDSNVKLEQLWKTNKQERKKINNKFCEQYITTFQKFDMDVQKFNEEQEKSNNYQKEQQALK
 LSKCSQSQTLEAIKDMHENYMEGLMNLNTNNYNMLFDVDGELRKEMSVFKKDLMKHTLKYSSFPSSD


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI



[View online »](#)

Cloning Scheme:



ACCN: NM_011725

ORF Size: 624 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_011725.4](#)

RefSeq Size: 1056 bp

RefSeq ORF: 627 bp


Locus ID: 22441

UniProt ID: [P05531](#)

Cytogenetics: X A5

MW: 24.6 kDa

Product images:



Circular map for MR202215