

Product datasheet for **MR20225A1V**

Mouse Ric3 (BC059258) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Ric3 (BC059258) AAV Particle
Tag: Myc-DDK
Symbol: Ric3
Synonyms: E130307J04Rik
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202225 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAAAAATTAATCAACAGAGTTGGACCTAATGGTGAGAGCAGAGCACAGGCTGTGACTTCTGACCAAG
AGAAACGATTACTGCATCAGCTCCGAGAAATCACCAGGGTCATGAAAGAAGGCAAGTTCATCGACACATC
TCCAGAGAAGGAAGCTGAGGAAGCCCCATACATGGAGGACTGGGAAGGTTACCCAGAAGAGACTTATCCA
ATATATGACCTTTCGGATGGTATCAAGCGTAGGCAAGAAACAATCCTGGTGGATTACCCTGACCTAAAAG
AGCCCTCTGCTGAAGAAATTGCTGAACAAATGGGAGAGATAGAAGAGGAAGGATCAGAGCGTTTGAGCTG
GGATCATTTGCCCACTGACCCTGGAGCCAGAAAGATAATTCTGTTGCCCGTGTGATCCCAAACCAGAA
TCATGCTCCTGCTGTGTTTCATGAAGAGGAGGATCCTGCTGTCTTGGCAGAGAATGCTGGTTTCAGTGCAG
ATGGCTACAGTGAGCAAGAGGAAGCCACCAAGGAAAACCTGCCCCAGGACTTTACAAATGAAGGGCTGGG
AGTCAGCACTGATAATGCACACGTGGTGGTATGTTGAGGAAGCGCAACCCCCAGGGTTTTGAG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202225 protein sequence
Red=Cloning site Green=Tags(s)

MEKLINRVGPNGESRAQAVTSDQEKRLHLQLREITRMKEGKFIDTSPEKEAEEAPYMEDWEGYPEETYP
IYDLSDGIKRRQETILVDYDLKEPSAEEIAEQMGEIEEEGSERLSWDHLPTDPAQKDNSVAPCDPKPE
SCSCCVHEEDPAVLAENAGFSADGYSEQEEATKENLPQDFTNEGLGVSTDNAHVGGMLRKNPQGFE

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Mouse



[View online »](#)

Serotype:	AAV-2
ACCN:	BC059258
ORF Size:	624 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC059258 , AAH59258
RefSeq Size:	5045 bp
RefSeq ORF:	626 bp
Locus ID:	320360
Cytogenetics:	7 E3
MW:	23.3 kDa