

## Product datasheet for **MR202228L3V**

### Gltp (NM\_019821) Mouse Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gltp (NM_019821) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Gltp
Synonyms:	1110001F24Rik; C76925; C77564
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_019821
ORF Size:	630 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202228).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_019821.2</a> , <a href="#">NP_062795.2</a>
RefSeq Size:	1622 bp
RefSeq ORF:	630 bp
Locus ID:	56356
UniProt ID:	<a href="#">Q9JL62</a>
Cytogenetics:	5 F



[View online »](#)

**Gene Summary:**

Accelerates the intermembrane transfer of various glycolipids. Catalyzes the transfer of various glycosphingolipids between membranes but does not catalyze the transfer of phospholipids. May be involved in the intracellular translocation of glucosylceramides (By similarity).[UniProtKB/Swiss-Prot Function]