

Product datasheet for **MR202232L3V**

Pgpep1 (NM_023217) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pgpep1 (NM_023217) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pgpep1
Synonyms:	2810003H13Rik; Gm19458; PAP-I; Pcp; PGP; PGP-I; PGPI
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_023217
ORF Size:	630 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202232).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_023217.2
RefSeq Size:	4882 bp
RefSeq ORF:	630 bp
Locus ID:	66522
UniProt ID:	Q9ESW8
Cytogenetics:	8 B3.3
Gene Summary:	Removes 5-oxoproline from various penultimate amino acid residues except L-proline. [UniProtKB/Swiss-Prot Function]



[View online »](#)