

Product datasheet for MR202247

Mrpl21 (BC028768) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mrpl21 (BC028768) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Mrpl21
Synonyms: MGC41352
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR202247 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTGCTGCCATAGCGGCCTCGGCATTGCCGGGCGCCTTCGGGCGGCTGGTGTCTGTGTAGCCGTA
 GCATCTTGGCATCCAGGGATCCGGATCAGCCTCCCTCTGGTCTGCTTCTCGAAGTTCAATTCACAGAG
 CGCTTCATATCCACAAGGATATGTTCTAAGACCTCCCTGAGCTCACCACCCTGGCAAGAAGTGGTTCTG
 CCGGACCCTGTTGAAGAGACGAGACACCATGCGGAGGTTGTGAAGAGGGTGAATGAGCTGATCGCCACGG
 GCCAGTATGGCCGGCTCTTTGCTGTCGTGCACTTTGCCAGCCACCAATGGAAGGTGACCGCAGAGGACCT
 GATTCTAATTGAAAATGAGTTGGATATCAAGTGTGGAGAGAGGATTCGGCTGGAGAAGGTCCTGCTGGTT
 GGGCAGACAACCTCACACTCCTGGCAAGCCACTCCTCCGAAAGGAGCTTGACGAGTTGAAGCCACAG
 TCATTGAAAAGACAGAATCATGGCCAAAATCAACATGAAATTTAGGAAAAGGAAAACTTCAGGAAGAA
 AAAAATCATTGTGAACCCTCAGACCATCCTCCGATTAACACCATTGAGATTGCACCTCGTTTGCTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202247 protein sequence
 Red=Cloning site Green=Tags(s)

MAAIIAASALPGAFGRLLVSVCSRSILASQSGSASLWSASRRFNSQSASYPQGYVPKTSLSPPWQEVVL
 PDPVEETRHHAEVVKRVNELIATGQYGRLLFAVVHFASHQWKVTAEDLIL IENELDIKCGERIRLEKVLV
 GADNFTLLGKPLLRKELVRVEATVIEKTESWPKINMKFRKRKNFRKKKIIVNPQILRINTIEIAPRL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: BC028768

ORF Size: 627 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [BC028768](#), [AAH28768](#)

RefSeq Size: 804 bp


RefSeq ORF: 629 bp

Locus ID: 353242

Cytogenetics: 19 A

MW: 23.4 kDa

Product images:



Circular map for MR202247