

Product datasheet for **MR202248**

Clec4n (NM_020001) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Clec4n (NM_020001) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Clec4n
Synonyms: Clec6a; Clecsf10; Nkcl
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR202248 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTGCAGGAAAGACAATCCCAAGGGAAGGGAGTCTGCTGGACCCTGAGACTCTGGTCAGCTGCTGTGA
 TTTCCATGTTACTCTTGAGTACCTGTTTCATTGCGAGCTGTGTGGTGACTTACCAATTTATTATGGACCA
 GCCAGTAGAAGACTATATGAACCTCACACATACCATTCCAGTCTCACCTGCTTCAGTGAAGGGACTATG
 GTGTCAGAAAAAATGTGGGGATGCTGCCAAATCACTGGAAGTCATTTGGCTCCAGCTGCTACCTCATTT
 CTACCAAGGAGAACTTCTGGAGCACCAGTGAGCAGAACTGTGTTTCAGATGGGGGCTCATCTGGTGGTGAT
 CAATACTGAAGCGGAGCAGAAATTCATCACCCAGCAGCTGAATGAGTCACTTTCTTACTTCTGGGTCTT
 TCGGATCCACAAGGTAATGGCAAATGGCAATGGATCGATGATACTCCTTTTCAGTCAAAATGTCAGTTCT
 GGCACCCCATGAACCAATCTCCAGAAGAGCGGTGTGTTTCAATAGTTTACTGGAATCCTTCGAAATG
 GGGCTGGAATGATGTTTTCTGTGATAGTAAACACAATTCAATATGTGAAATGAAGAAGATTTACCTA

A**CGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202248 protein sequence
 Red=Cloning site Green=Tags(s)


MVQERQSQGKGVCTRLWSAAVISMLLLSTCFIASCVVTYQFIMDQPSRRLYELHTYHSSLTCFSEGT
 VSEKMWGCCPNHWKSGFGSSCYLISTKENFWSTSEQNCVQMG AHLVVINTEAEQNFITQQLNESLSYFLGL
 SDPQNGKQWIDDPFSQNVRFWHPHEPNLPEERCVSIVYWNPSKWGWNDVFCDSKHNSICEMKKIYL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Restriction Sites: Sgfl-MluI

Cloning Scheme:



ACCN: NM_020001

ORF Size: 627 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_020001.2](#)

RefSeq Size: 1245 bp

RefSeq ORF: 630 bp

Locus ID: 56620


UniProt ID: Q9JKF4

Cytogenetics: 6 58.3 cM

MW: 24.3 kDa

Gene Summary: Binds high-mannose carbohydrates in a Ca(2+)-dependent manner. Functional receptor for alpha-mannans on *C.albicans* hyphae. Plays an important role in the host defense against *C.albicans* infection by inducing TH17 cell differentiation. Recognizes also, in a mannose-dependent manner, allergens from house dust mite and fungi, by promoting cysteinyl leukotriene production. Recognizes soluble elements from the eggs of *Shistosoma mansoni* altering adaptive immune responses. Transduces signals through an Fc receptor gamma chain /FCER1G and Syk-CARD9-NF-kappa-B-dependent pathway.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR202248