

## Product datasheet for **MR202248A1V**

### Mouse Clec4n (NM\_020001) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Mouse Clec4n (NM\_020001) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** Clec4n  
**Synonyms:** Clec6a; Clecsf10; Nkcl  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >MR202248 ORF sequence  
**Red**=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGTGCAGGAAAGACAATCCCAAGGGAAGGGAGTCTGCTGGACCCTGAGACTCTGGTCAGCTGCTGTGA  
TTTCCATGTTACTCTTGAGTACCTGTTTCATTGCGAGCTGTGTGGTGACTTACCAATTTATTATGGACCA  
GCCAGTAGAAGACTATATGAATTCACACATACCATTCAGTCTCACCTGCTTCAGTGAAGGGACTATG  
GTGTCAGAAAAATGTGGGATGCTGCCAAATCACTGGAAGTCATTTGGCTCCAGCTGCTACCTCATTT  
CTACCAAGGAGAACTTCTGGAGCACCAGTGAGCAGAACTGTGTTTCAGATGGGGGCTCATCTGGTGGTGAT  
CAATACTGAAGCGGAGCAGAATTTTCATCACCAGCAGCTGAATGAGTCATTTCTTACTTCTCTGGGTCTT  
TCGGATCCACAAGGTAATGGCAATGGCAATGGATCGATGATACTCCTTTCAGTCAAAATGTCAGGTTCT  
GGCACCCCATGAACCAATCTTCCAGAAGAGCGGTGTGTTTCAATAGTTTACTGGAATCCTTCGAAATG  
GGGCTGGAATGATGTTTTCTGTGATAGTAAACACAATTCAATATGTGAAATGAAGAAGATTTACCTA

**ACGCGT**ACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR202248 protein sequence  
**Red**=Cloning site **Green**=Tags(s)

MVQERQSQGKGVCTLRLLWSAAVISMLLLSTCFIASCVVTYQFIMDQPSRRLYELHTYHSSLTCFSEGT  
VSEKMWGCCPNHWKSFSGSSCYLISTKENFWSTSEQNCVQMG AHLVVINTEAEQNFITQQLNESLSYFLGL  
SDPQNGKQWIDDPFSQNVRFWHPHEPNLPEERCVSIVYWNPSKWGWNDVFCDSKHNSICEMKKIYL

**TR**TRPLE**QKL**ISEEDLAANDILDYKDDDDKV

**Species:** Mouse



[View online »](#)

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Serotype:</b>     | AAV-2                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>ACCN:</b>         | NM_020001                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>ORF Size:</b>     | 627 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Buffer:</b>       | PBS with 0.001% Pluronic F68                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Stability:</b>    | AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus. |
| <b>RefSeq:</b>       | <u><a href="#">NM_020001.1</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>RefSeq Size:</b>  | 1245 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>RefSeq ORF:</b>   | 630 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Locus ID:</b>     | 56620                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>UniProt ID:</b>   | <u><a href="#">Q9JKF4</a></u>                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Cytogenetics:</b> | 6 58.3 cM                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>MW:</b>           | 24.3 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |