

Product datasheet for **MR202264L4V**

Fgf21 (NM_020013) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgf21 (NM_020013) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fgf21
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_020013
ORF Size:	630 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202264).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_020013.4 , NP_064397.1
RefSeq Size:	947 bp
RefSeq ORF:	633 bp
Locus ID:	56636
UniProt ID:	Q9JIN1
Cytogenetics:	7 B3
Gene Summary:	Stimulates glucose uptake in differentiated adipocytes via the induction of glucose transporter SLC2A1/GLUT1 expression (but not SLC2A4/GLUT4 expression). Activity probably requires the presence of KLB.[UniProtKB/Swiss-Prot Function]



[View online »](#)