

Product datasheet for **MR202286L4V**

Fbxo44 (BC028884) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fbxo44 (BC028884) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Fbxo44
Synonyms:	Fbx6a, FBG3, FBX30, 5730411K09
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC028884
ORF Size:	633 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202286).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC028884 , AAH28884
RefSeq Size:	974 bp
RefSeq ORF:	635 bp
Locus ID:	230903
Cytogenetics:	4 E2
Gene Summary:	Substrate-recognition component of the SCF (SKP1-CUL1-F-box protein)-type E3 ubiquitin ligase complex.[UniProtKB/Swiss-Prot Function]



[View online »](#)