

## Product datasheet for **MR202305L3V**

### **Ndufs8 (NM\_144870) Mouse Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Ndufs8 (NM_144870) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Ndufs8
Synonyms:	BC021616; CI-23kD; TYKY
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144870
ORF Size:	639 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202305).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_144870.3</a> , <a href="#">NP_659119.2</a>
RefSeq Size:	1057 bp
RefSeq ORF:	639 bp
Locus ID:	225887
UniProt ID:	<a href="#">Q8K3J1</a>
Cytogenetics:	19 A



[View online »](#)

**Gene Summary:**

Core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I) that is believed to belong to the minimal assembly required for catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone (By similarity). [UniProtKB/Swiss-Prot Function]