

Product datasheet for **MR202306L1V**

Pdcd10 (NM_019745) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pdcd10 (NM_019745) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Pdcd10
Synonyms:	2410003B13Rik; Ccm3; Tfa15; Tfar15
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_019745
ORF Size:	636 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202306).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_019745.2 , NP_062719.2
RefSeq Size:	1923 bp
RefSeq ORF:	639 bp
Locus ID:	56426
UniProt ID:	Q8VE70
Cytogenetics:	3 E3



[View online »](#)

Gene Summary:

Promotes cell proliferation. Modulates apoptotic pathways. Increases mitogen-activated protein kinase activity and STK26 activity. Important for cell migration, and for normal structure and assembly of the Golgi complex (By similarity). Important for KDR/VEGFR2 signaling. Increases the stability of KDR/VEGFR2 and prevents its breakdown. Required for normal cardiovascular development. Required for normal angiogenesis, vasculogenesis and hematopoiesis during embryonic development (By similarity).[UniProtKB/Swiss-Prot Function]