

Product datasheet for **MR202322A1V**

Mouse Nabp2 (NM_027257) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Nabp2 (NM_027257) AAV Particle
Tag: Myc-DDK
Symbol: Nabp2
Synonyms: 2610036N15Rik; Obfc2b; SSB1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202322 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGACGACCGAGACCTTCGTGAAGGATATCAAGCCCGGGCTCAAGAATCTGAATCTTATCTTCATTGTAC
 TGGAGACAGGCCGGGTGACCAAGACAAAGGACGGACACGAGGTTTCGGACCTGTAAAGTGGCAGACAAAAC
 GGGCAGCATCAACATCTCAGTATGGGACGATGTGGCAACCTGATCCAACCTGGGGACATTATCCGACTC
 ACCAAAGGGTATGCTTCAGTGTTCAAAGTTGTCTGACACTGTACACTGGCCGTGGGGCGATCTGCAGA
 AGATTGGAGAATTCTGTATGGTTTATTCGAAGTCCCTAATTTTCAGTGAGCCCAATCCAGAGTATAACAC
 GCAGCAGGCGCCCAACAAATCGGTGCAGAACAAATGACAACAGCCCTACTGCCCCACAAGCGACACGCGGA
 CCCCCTGCTGCCTCTCCAGCTTCGAGAACGAGAACGAAATGGACTGAGCACCCAGCTAGGTCTGTGG
 GTGGCCCCCATCCTTCTCACACCCCTCACATCCCCCAGCACCCGAATCACTCGAAGTCAACCAACCA
 CACGCCTTCGGCCCTCTGGCCCTCCAGCAACCCTGTCAGCAACGGCAAGAAACCCGAAGAAGCAGC
 AAGAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202322 protein sequence
 Red=Cloning site Green=Tags(s)

MTTETVFKDIKPLKLNLIIFIVLETGRVTKTKDGHEVRTCKVADKTGSINISVWDDVGNLIQPGDIIRL
 TKGYASVFKGCLTLYTGRGDLQKIGEFCEVYSEVPNFSEPNPEYNTQQAPNKSQNNDSPTAPQATTG
 PPAASPAENQNGNLSTQLGPVGGPHPSHTPSHPSTRITRQPNHTPSGPPGPSSNPVSNGKETRRSS
 KR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_027257
ORF Size:	636 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_027257.1
RefSeq Size:	1106 bp
RefSeq ORF:	639 bp
Locus ID:	69917
UniProt ID:	Q8R2Y9
Cytogenetics:	10 D3
MW:	22.6 kDa