

Product datasheet for **MR202344A1V**

Mouse Tnfrsf19 (BC030062) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Tnfrsf19 (BC030062) AAV Particle
Tag: Myc-DDK
Symbol: Tnfrsf19
Synonyms: TAJ, Troy, TAJ-ALPHA
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202344 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCACTCAAGGTCCTACCTCTACACAGGACGGTCTTCGCTGCCATTCTCTTCTACTCCACCTGG
CATGTAAGTGAGTTGCGAAACCGGAGATTGCAGGCAGCAGGAATTCAGGATCGATCTGGAACTGTGT
CCTCTGCAACAGTGCGGACCTGGCATGGAGTTGTCCAAGGAATGTGGCTTCGGCTATGGGGAGGATGCA
CAGTGTGTGCCCTGCAGGCCGCACCGGTTCAAGGAAGACTGGGGTTTCCAGAAGTGAAGCCATGTGCGG
ACTGTGCGCTGGTGAACCGCTTTCAGAGGGCCAACTGCTCACACACCACTGATGCTGTCTGCGGGGACTG
CCTGCCAGGATTTTACCGGAAGACCAAACTGGTTGTTTCAAGACATGGAGTGTGTGCCCTGCGGAGAC
CCACCTCCTCCCTACGAACCACACTGTACCAGCAAGGTGAACCTTGTGAAGATCTCCTCCACCGTCTCCA
GCCCTCGGGACACGGCGCTGGCTGCCGTATCTGCAGTGTCTGGCCACGGTGTCTGCGCCCTGCTCAT
CCTGTGTGTCTACTGCAAGAGGCAGTTTCATGGAGAAGAAACCCAGCTGTAAGCTCCCATCCCTCTGT
CTCACTGTGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202344 protein sequence
Red=Cloning site Green=Tags(s)

MALKVLPLHRTVLF AAILFLLHLACKVSCETGDCRQQEFKDRSGNCVLCKQCGPGMELSKCEGFGYGEDA
QCVPCRPHRFKEDWGFQKCKPCADCALVNRFRANCSTSDAVCGDCLPGFYRKTLVGFQDMECVPCGD
PPPPYEPHCTSKVNLVKISSTVSSPRDTALAAVICSALATVLLALLILCVIYCKRQFMKKPCKLPSLC
LTVK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	BC030062
ORF Size:	642 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC030062 , AAH30062
RefSeq Size:	1271 bp
RefSeq ORF:	644 bp
Locus ID:	29820
Cytogenetics:	14 D1
MW:	23.6 kDa