

Product datasheet for **MR202344L3V**

Tnfrsf19 (BC030062) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tnfrsf19 (BC030062) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Tnfrsf19
Synonyms:	TAJ, Troy, TAJ-ALPHA
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	BC030062
ORF Size:	642 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202344).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC030062 , AAH30062
RefSeq Size:	1271 bp
RefSeq ORF:	644 bp
Locus ID:	29820
Cytogenetics:	14 D1
Gene Summary:	Can mediate activation of c-Jun and NF-kappa-B. May promote caspase-independent cell death (By similarity). Isoform 2 and isoform 3 may act as decoy receptors.[UniProtKB/Swiss-Prot Function]



[View online »](#)