

Product datasheet for **MR202371A1V**

Mouse Adamts15 (BC034843) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Adamts15 (BC034843) AAV Particle
Tag: Myc-DDK
Symbol: Adamts15
Synonyms: MGC41037
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202371 representing BC034843
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCAATTCACCCTGACCCTGCTCTTCACTGCCTCTTGGTAACCCTGTGCTCCCTCAAACCCCGCCTCT
TCCCGCTTCTCCGTTCTTGCTGTGGACTCTCTGAGCTGTGATTGAGGGGTGCTGCCAGGGTCCGGG
AGAGTGGACGCCCTGGGGTCTTGGAGCAGGTGTTCCAGCTCTTGTGGTCGGGGCGTCTCAGTGCCGAGA
CGCCAATGTGTCCGACTTCCCGAGAGGAACCTTGTGGGGCGACAGCCATGAATACCGGGTGTGCCAGC
TGCCGGACTGCCCTCCAGGGACACACCATTCGAGACCTGCAGTGTTCTCTACAATGCCACCTGT
CCTGGACACCCAGAAAACCTACCACTGGGTGCCCTTCCACGGGGCGCCCAATGTGTGTGACCTCAACTGC
CTAGCTGAAGGCCACGCCTTCTACCACAGCTTTGGCCGCGTGCTGGATGGCACCCCTGCACCCCGGGTA
CCCAGGGTCTCTGCGTGGCAGGCCGCTGTCTGAGCGCAGGCTGTGACGGGATACTGGGCTCCGGCACCT
CGAGGACCGCTGCGGCCCTCTGTGGTGGCGCTAACGACTTTGCCTGTTCTGTCAGCGCGTGTCCGCGAC
GCCGCTATCGGTCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202371 representing BC034843
Red=Cloning site Green=Tags(s)

MAIHPDPALHCLLVTLCSLKPRLLPFLWTLSCDLRGAAQGPGEWTPWGSWSRCSSSCGRGVSRR
RQCVRLPGEELCWGDSHEYRVCLPDCPPGTTTFRDLQCSLYNGHPVLDTKTYQWVPFHGAPNVCDLNC
LAEGHAFYHSFGRVLDGTPGTQGLCVAGRCLSGCDGILSGTLEDRCGLCGGANDSCLFVQRVFRD
AAIGP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	BC034843
ORF Size:	645 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	BC034843
RefSeq Size:	2099 bp
RefSeq ORF:	647 bp
Locus ID:	66548
Cytogenetics:	10 C1
MW:	76.9 kDa