

Product datasheet for **MR202381L3V**

Sspn (NM_010656) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sspn (NM_010656) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Sspn
Synonyms:	Krag
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_010656
ORF Size:	648 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202381).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_010656.2 , NP_034786.1
RefSeq Size:	4408 bp
RefSeq ORF:	651 bp
Locus ID:	16651
UniProt ID:	Q62147
Cytogenetics:	6 77.7 cM



[View online »](#)

Gene Summary:

Component of the dystrophin-glycoprotein complex (DGC), a complex that spans the muscle plasma membrane and forms a link between the F-actin cytoskeleton and the extracellular matrix. Preferentially associates with the sarcoglycan subcomplex of the DGC (By similarity). [UniProtKB/Swiss-Prot Function]