

Product datasheet for **MR202388L4V**

Mob3b (NM_178061) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Mob3b (NM_178061) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Mob3b
Synonyms:	8430436F23Rik; A430018A01Rik; Mobkl2b
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_178061
ORF Size:	651 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202388).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_178061.4
RefSeq Size:	6030 bp
RefSeq ORF:	651 bp
Locus ID:	214944
UniProt ID:	Q8VE04
Cytogenetics:	4 16.86 cM
Gene Summary:	Modulates LATS1 expression in the Hippo signaling pathway which plays a pivotal role in organ size control and tumor suppression by restricting proliferation and promoting apoptosis.[UniProtKB/Swiss-Prot Function]



[View online »](#)