

Product datasheet for **MR202397L4V**

Rassf4 (BC023245) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rassf4 (BC023245) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rassf4
Synonyms:	3830411C14Rik; AD037; AI586137
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	BC023245
ORF Size:	648 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202397).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	BC023245 , AAH23245
RefSeq Size:	2817 bp
RefSeq ORF:	650 bp
Locus ID:	213391
Cytogenetics:	6 E3
Gene Summary:	Potential tumor suppressor. May act as a KRAS effector protein. May promote apoptosis and cell cycle arrest.[UniProtKB/Swiss-Prot Function]



[View online »](#)