

Product datasheet for **MR202445L3V**

Rab27b (NM_030554) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Rab27b (NM_030554) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Rab27b
Synonyms:	2310021G14Rik; B130064M09Rik
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_030554
ORF Size:	657 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202445).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030554.3 , NP_085031.3
RefSeq Size:	6977 bp
RefSeq ORF:	657 bp
Locus ID:	80718
UniProt ID:	Q99P58
Cytogenetics:	18 43.89 cM
Gene Summary:	Plays a role in NTRK2/TRKB axonal anterograde transport by facilitating the association of NTRK2/TRKB with KLC1. May be involved in targeting uroplakins to urothelial apical membranes.[UniProtKB/Swiss-Prot Function]



[View online »](#)