

Product datasheet for **MR202455A1V**

Mouse Polr3gl (NM_027241) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Polr3gl (NM_027241) AAV Particle
Tag: Myc-DDK
Symbol: Polr3gl
Synonyms: 2010003I19Rik
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202455 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCAGCCGGGGTGGGGGCCGGGGTCGTGGCCGGGGCCAGTTGACCTTCAACATGGAGGCTGTGGGCA
TTGGGAAGGGTGATGCTTTGCCCCACCCACTCTCCAGCCCTCTCCACTCTTCTCCTCCCTGGAGTTCCA
CCCAGTGCCTCTGCCGCGAGGAGGAGGGGGAGTATGTCCTGGCACTGAAGCAAGAGCTGCGCGGGGCC
ATGAGGCAGCTCCCCTACTTCATCCGGCCAGCCGTCCCCAAGAGAGATGTGGAACGTTACTCAGACAAGT
ATCAGATGTCTGGGCTATTGACAACGCCATCGATTGGAACCCGTGATTGGCGGCGACTCCCCAGTGAGCT
CAAGATTCGAGTGCGGAAAGTACAGAAGGAGCGGACCACCATTATCCTTCCCAAGAGGCCCCCTAAGAGC
ACAGACGATAAGGAGGAGACGATACAGAACTAGAGACGCTAGAGAAGAAGGAGGAGGAAGTGACTTCGG
AGGAGGATGAAGAGAAAGAAGAAGAGGAGGAGAAGGAAGAGGGGAGGAAGAGGAATACGATGAAGAGGA
GCATGAGGAGGAACTGATTACATCATGTATTTTGACAATGGAGAGGACTTTGGAGGGGACAGCGAT
GACAATATGGATGAGGCCATATACTGAAAACTGGGCATAGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202455 protein sequence
Red=Cloning site Green=Tags(s)

MASRGGGRGRGRGQLTFNMEAVGIGKGDALPPPTLQPSPLFPPLEFHPVPLPAGEEGEYVLALKQELRGA
MRQLPYFIRPAVPKRDVERYSDKYQMSGPIDNAIDWNPDRRLPSELKIRVRKVQKERTTIILPKRPPKS
TDDKEETIQLETLEKKEEEVTSSEDEEKEEEEEKEEEEEYDEEEHEEETDYIMSYFDNGEDFGGDS
DNMDEAHILKNLGID

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_027241
ORF Size:	675 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_027241.1</u>
RefSeq Size:	1427 bp
RefSeq ORF:	657 bp
Locus ID:	69870
UniProt ID:	<u>Q8R0C0</u>
Cytogenetics:	3 F2.1
MW:	25.9 kDa