

Product datasheet for **MR202462**

Mageh1 (NM_023788) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mageh1 (NM_023788) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Mageh1
Synonyms:	2010107K23Rik; APR-1; APR1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202462 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCCTCGGGGACGGAAGAGTCGGCGCCGCCGGAACGCAAAGGCAGCTGAAGAGAATCGCAACAATCGCA
 AGAGCCAGGCCTCAGAGGCTTCAGAGACCCGATGGCGGCTTCTGTAGCCCCGAGCACACCCGAAGAATA
 CCTGAGCGGCCCGGAGGAAGACACAAGCACCTGGAGAAGGCCTCCAGTACCCCTTCAGAAGCTTCGAGC
 ACTGCCCTAGTGAAAAGCCGGTTACCCGAGCAATTTTCAGGGCACCAAGAAAAGTCTCCTAATGTCCA
 TATTAGCCCTCATCTTCATCATGGGCAACAGCGCAAAGGAGGCCCTGGTGTGAAAGTGCTGGGGAAGTT
 AGGGATGCAGCCTGGGAGGCAGCACAGCATCTTTGGAGATCCGAAGAAGGTCGTTACAGAAGAGTTTGTG
 CGCAGAGGGTATCTGATTATAAGCCAGTGCCCGCAGCAGTCCAGTGGAGTATGAGTTCTCTGGGGCC
 CTCGAGCACACGTGGAATCTAGCAAGCTGAAAGTCATGCATTTTGTGGCAAGAGTTCGGAACCGATGCTC
 CAAAGACTGGCCATGTAATATGACTGGGATTCAGATGATGATGCAGAGGTTGAGGCTATTCTCAATTCC
 GGTGCTAGGGTTATTCTGCTCCA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >MR202462 protein sequence
 Red=Cloning site Green=Tags(s)



MPRGRKSRRRRNAKAAEENRNRKSQASESETPMAASVAPSTPEEYLSGPEEDTSTLEKASSTPSEASS
 TALVQKPVTRSNFQGTKKSLMSILALIFIMGNSAKEALVVKVLGKLGMPGRQHSIFGDPKKVVTEEFV
 RRGYLIYKPVPRSSPVEYEFWGPRAHVESSKLKVMHFVARVRNRCSDKDWPCNYDWDSDDDAEVEAILNS
 GARGYSAP

TRTRPLEOKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM 023788

ORF Size: 654 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_023788.3](#), [NP_076277.1](#)

RefSeq Size: 1398 bp

RefSeq ORF: 657 bp


Locus ID: 75625

UniProt ID: [Q9NWX9](#)

Cytogenetics: X F3

MW: 24.3 kDa

Product images:



Circular map for MR202462