

## Product datasheet for MR202464L4

### Taf10 (NM\_020024) Mouse Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Taf10 (NM_020024) Mouse Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Taf10
Synonyms:	30kDa; AU041226; Taf2h; TAFII30
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202464).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 24075

UniProt ID: [Q8K0H5](#)

Cytogenetics: 7 E3

**Gene Summary:** TAFs are components of the transcription factor IID (TFIID) complex, PCAF histone acetylase complex and TBP-free TAFII complex (TFTC). TFIID is a multimeric protein complex that plays a central role in mediating promoter responses to various activators and repressors. May regulate cyclin E expression.[UniProtKB/Swiss-Prot Function]