

Product datasheet for **MR202479L4V**

Trappc4 (NM_021789) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Trappc4 (NM_021789) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Trappc4
Synonyms:	1500017G03Rik; AI303617; HSPC172; PTD009; Sbd; Sbdn; Sdcbp2; TRS23
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_021789
ORF Size:	660 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202479).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021789.2 , NP_068561.1
RefSeq Size:	1196 bp
RefSeq ORF:	660 bp
Locus ID:	60409
UniProt ID:	Q9ES56
Cytogenetics:	9 A5.2
Gene Summary:	May play a role in vesicular transport from endoplasmic reticulum to Golgi. May play a role in dendrite postsynaptic membrane trafficking.[UniProtKB/Swiss-Prot Function]



[View online »](#)