

Product datasheet for **MR202482A1V**

Mouse Cideb (NM_009894) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Cideb (NM_009894) AAV Particle
Tag: Myc-DDK
Symbol: Cideb
Synonyms: 1110030C18Rik; AI790179; CIDE-B
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202482 ORF sequence
Red=Cloning site **Blue**=ORF **Green**=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGTACCTTTCAGCCTTCAACCCCAATGGCCTGCTAAGGTCAGTATCCACTGTGAGCTCGGAGTTAA
GCCGTAGGGTCTGGAAGCTCAGCTCCTCCACCTCAGCGACCTTCCGTGTCTGTGATCATAAGCGGACAGT
CCGGAAGGACTGACAGCTGCCAGCCTCCAAGAACTGCTAGATAAGGTCCTGGAGACCTTGCTGCTACGT
GGAGTGCTAACACTGGTCCTGGAGGAGGATGGGACTGCTGTGGACAGTGAGGACTTCTCCAGCTGCTGG
AGGACGACACGTGCTTGATGGTGCTTGAGCAGGGCCAGAGCTGGAGCCCCAAGAGTGGGATGTTGCATA
CGGCCCTAGGACGGGAGAAGCCAAAACACAGCAAGGACATCGCCCGCATCACCTTCGATGTGTACAAGCAA
AATCCCCGAGACCTCTTTGGCAGCCTCAACGTGAAAGCTACATTCTATGGACTCTACTCCATGAGCTGTG
ATTTCCAAGGAGTTGGCCCTAAAAGAGTACTCAGGGAGCTCCTCCGTTGGACTTCCTCGCTGCTGCAAGG
CCTGGGCCATATGCTGCTGGGCATCTCTCCACCCTTCGCCATGTGGTGGAGGGGGCTGATCGATGGCAG
TGGCACGGGCAGAGACACCTCCACTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202482 protein sequence
Red=Cloning site **Green**=Tags(s)

MEYLSAFNPNGLLRSVSTVSSELSRRVWNSAPPPQRPFRVCDHKRTVRKGLTAASLQELLDKVLLETLRLR
GVLTLVLEEDGTAVDSEDFQLLEDDTCLMVLEQQQSWSPKSGMLSYGLGREPKHKSIDIARITFDVYKQ
NPRDLFGSLNVKATFYGLYSMSCDFQGVGPKRVLRELLRWTSSLLQGLGHMLLGISSTLRHVVEGADRWQ
WHGQRHLHS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_009894
ORF Size:	660 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_009894.1
RefSeq Size:	1186 bp
RefSeq ORF:	660 bp
Locus ID:	12684
UniProt ID:	O70303
Cytogenetics:	14 C3
MW:	24.8 kDa