

Product datasheet for **MR202486L3V**

Cldn6 (NM_018777) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cldn6 (NM_018777) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Cldn6
Synonyms:	AL024037
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018777
ORF Size:	660 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202486).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018777.2
RefSeq Size:	1522 bp
RefSeq ORF:	660 bp
Locus ID:	54419
UniProt ID:	Q9Z262
Cytogenetics:	17 A3.3



[View online »](#)

Gene Summary:

This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. The protein encoded by this gene is essential for blastocyst formation in preimplantation mouse embryos, and is involved in and is crucial for the formation and maintenance of the epidermal permeability barrier. This gene is adjacent to another family member Cldn9 on chromosome 17. [provided by RefSeq, Aug 2010]