

Product datasheet for **MR202502L3V**

Creg1 (NM_011804) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Creg1 (NM_011804) Mouse Tagged ORF Clone Lentiviral Particle
Symbol:	Creg1
Synonyms:	AA755314; Creg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_011804
ORF Size:	663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(MR202502).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_011804.2
RefSeq Size:	2138 bp
RefSeq ORF:	663 bp
Locus ID:	433375
UniProt ID:	O88668
Cytogenetics:	1 H2.3



[View online »](#)

Gene Summary:

May contribute to the transcriptional control of cell growth and differentiation. Antagonizes transcriptional activation and cellular transformation by the adenovirus E1A protein. The transcriptional control activity of cell growth requires interaction with IGF2R (By similarity). [UniProtKB/Swiss-Prot Function]