

Product datasheet for **MR202513A1V**

Mouse U2af1l4 (NM_170760) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse U2af1l4 (NM_170760) AAV Particle
Tag: Myc-DDK
Symbol: U2af1l4
Synonyms: AA407033; AF419339; AI451269; AW553050; U2af26
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202513 representing NM_170760
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCTGAATATTTAGCTTCGATATTCGGGACTGAGAAGGACAAGGTTAACTGCTCTTTTACTTTAAGA
TTGGAGCCTGCCGGCACGGGGACCGGTGCTCCCGACTTCACAACAAACCGACTTTCAGCCAGACCATAGT
CCTGCTCAACTTGTACCGGAATCCACAGAACACAGCCAACTGCAGATGGATCACACTGTCACGTGAGC
GACGTGGAGGTGCAAGAACACTATGATAACTTCTTTGAGGAGGTATTCACAGAACTGCAGGAGAAGTATG
GAGGATTGAAGAGATGAATGTGTGTGACAACCTCGGGGACCACCTCGTGGGCAATGTCTACGTTAAGTT
CGGGCGGAGGAGGATGCAGAGCGGGCTGTAGCGGAACCTCAATAACCGCTGGTTCAACGGGCAGGCTGTG
CATGCCGAGCTGTCTCCTGTCACTGACTTCCGAGAGTCCTGCTGCCGGCAGTATGAGATGGGGGAATGCA
CCCGAGGTGGCTTCTGCAACTTTATGCACCTACGGCCCATATCTCGGAACTTGCGCCGGCAGCTCTATGG
GCGAGGACCCAGGCATAGGTACCTCCAAGGTCCACACAGGTACCGTCCCCGAGAAAGGAACCGACGT
CGTCCCCAGACCACCGCATGGTCGCTTC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202513 representing NM_170760
Red=Cloning site Green=Tags(s)

MAEYLASIFGTEKDKVNCSFYFKIGACRHGDRCSRLHNKPTFSQTIVLLNLYRNPQNTAQTAGDSHCHVS
DVEVQEHYDNFFEEVFTLQEKYGEIEEMNVCDNLGDHLVGNVYVKFRREEDAERAVAELENNRWFNGQAV
HAELSPVTDFFRESCCRQYEMGECTRGGFCNFMHLRPISRNLRRQLYGRGPRHRSPPRSHTGHRPRERNRR
RSPDHRHGRF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_170760
ORF Size:	660 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_170760.3 , NP_739566.1
RefSeq Size:	864 bp
RefSeq ORF:	663 bp
Locus ID:	233073
UniProt ID:	Q8BGJ9
Cytogenetics:	7 B1
MW:	25.8 kDa