

Product datasheet for MR202519

Lhfp1 (NM_178358) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Lhfp1 (NM_178358) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Lhfp1
Synonyms:	Lhfp
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR202519 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGGAACAGCCTTACTATGGTAGGCACCTTCTGGGCCCTCCTGTCCCTAGTTACCGCTGTGGCCAGTT
CTACCAGTTACTTCCTTCTTACTGGCTGTTGGATCACAGCTTGGAAAGCCAGTATCTTTCAGCACATT
CAGGAGATGCAACTACCCTGTGCGGGGAGATGGGCACAATCTGATCATGGTGAAGAATGTGGGCGCTAT
GCCAGCTTCAACCCATACCAAGCCTGGCCTGGCAGATGTGCACGGTGGTACAGGTGCTGGCTGTGCTC
TGTTGCTTCTGGTGGCATTGGCTGCTGTCTAGGATGCTGCATGGAAGAGCTCATCTCTAGAATGATGGG
GGTGGATGGGAGCTGCACAGTTTGTAGGAGGCTGCTGATAAGCGCAGGCTGTGCATTATATCCTTTA
GGATGGAACAGCCCGGAGGTAATGCAAACATGCGGGAATGTCTCCAATCAATTTTCAGTTAGGAACCTGTA
GACTCGGCTGGGCTATTACTGTGCTGGGGTGGAGCTGCTGCAGCCATGTTGATCTGCACCTGGCTCTC
TTGCTTTGCTGGAAGAAACCCCAAGCCAGTTATGTTGGTGGAGAACATCATGAGAAATACCAATTCTTAT
GCCTTGGAGCTTGATCACTGCCTCAAACCT

ACGCGTACGCGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR202519 protein sequence
 Red=Cloning site Green=Tags(s)

MRNSLTMVGTFWAFLSLVTAVASSTSYFLPYWLFQSGLGKPVSFSTFRRCNYPVRGDGHNLIMVEECGRY
 ASFTAIPSLAWQMCTVVTGAGCALLLLVALAAVLGCCMEELISRMMGRCMGAAQFVGLLISAGCALYPL
 GWNSPEVMQTCGNVSNQFQLGTCRLGWAYYCAGGAAAAMLICTWLSCFAGRNP KP VMLVENIMRNTNSY
 ALELDHCLKP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_178358

ORF Size: 663 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_178358.4](#)

RefSeq Size: 1827 bp

RefSeq ORF: 663 bp


Locus ID: 237091

UniProt ID: [Q80SV1](#)

Cytogenetics: X F2

MW: 23.8 kDa

Product images:



Circular map for MR202519