

Product datasheet for **MR202526A1V**

Mouse Tmem192 (NM_001163747) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Tmem192 (NM_001163747) AAV Particle
Tag: Myc-DDK
Symbol: Tmem192
Synonyms: 3110005G23Rik
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR202526 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGCGGCCGCGGCTGGAGGACAGTTCCTTGGATATCCTACAAAGTATGGATGACGACCCGCTTC
 TGGATACCCAGCCTCTCCACACCACTCTCTCAAGCTCAGTTTAGACCCAGATTCATCCTCTCCCCAC
 GGTCACTATAGCGAACCTTCTGTTGCTTATACATGTCGTGTTTGTGTTTAGCATTTTTAACGGGTGTG
 CCGTGTGTTGTATCCTAATCCAAGTGAAGACAAGTGCCAGAAAACCTACACGAGCCCGTTGAAGGTTTCA
 CGGCCATAATCCTTGGGAACTTATCCTGTGATTCTGCATTTGCTTTTGAACGCTACGTCCAGTATCA
 CCACAGGAAAGTCAGGAGCCGAGGCTATTCCAGATCTACCGCTCTACGAGGCATCTCAAAACACTGGCC
 CTCACCATCCACTCTTCTGTGAAGATTAGGAGGTTCAACAGAGCCAAGCCACTGCCTGATGTACTTGAAG
 AAGAGAAAATTTATGCTTACCCAGCAATACTGCCTCAGAGACTGGCTTCAGGACTGTTTCAAGCCTAGA
 AGAAATTGTTGAAAAGCAAGAGGACATAATTGTGTACCTGAAGCGACACAATGCTCTGCTGAGTAAGCGG
 TTGTTGGAGCTGGCAACTCAGCCAGCCAGGACT

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR202526 protein sequence
 Red=Cloning site Green=Tags(s)

MAAAGRLEDSSLDILQSMDDPLDQTPLPHSLQAQFRPRFHLPTVIANLLLLIHVVFVLAFLTG
 PCLYPNPTEDKCPENYTSPLKVQTAILGLILWLHLLFERYVQYHHRKVRSGYSQIYRSTRHLKTLA
 LTIHSSVKIRRFNRKPLPDVLEEEKIYAYPSNTASETGFRTVSSLEEIVEKQEDIIVYLKRHNALLSKR
 LLELATQPART

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


[View online »](#)

Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001163747
ORF Size:	663 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001163747.1</u> , <u>NP_001157219.1</u>
RefSeq Size:	1005 bp
RefSeq ORF:	666 bp
Locus ID:	73067
UniProt ID:	<u>Q9CXT7</u>
Cytogenetics:	8 B3.1
MW:	25.4 kDa