

Product datasheet for **MR202527L4V**

Trp53inp2 (NM_178111) Mouse Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Trp53inp2 (NM_178111) Mouse Tagged ORF Clone Lentiviral Particle |
| Symbol: | Trp53inp2 |
| Synonyms: | 1110029F20Rik; Tp53inp2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_178111 |
| ORF Size: | 666 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(MR202527). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_178111.3 , NP_835212.1 |
| RefSeq Size: | 3934 bp |
| RefSeq ORF: | 666 bp |
| Locus ID: | 68728 |
| UniProt ID: | Q8CFU8 |
| Cytogenetics: | 2 H1 |


[View online »](#)

Gene Summary:

Dual regulator of transcription and autophagy. Positively regulates autophagy and is required for autophagosome formation and processing. May act as a scaffold protein that recruits MAP1LC3A, GABARAP and GABARAPL2 and brings them to the autophagosome membrane by interacting with VMP1 where, in cooperation with the BECN1-PI3-kinase class III complex, they trigger autophagosome development. Acts as a transcriptional activator of THRA. [UniProtKB/Swiss-Prot Function]